● **TOYOBO BIOCHEMICAL PRODUCTS** ●
(Diagnostic Reagent Grade)

Blocking Peptide Fragment

from Microorganism

PREPARATION and SPECIFICATION

Appearance	: White amorphous powder, lyophilized
Grade	: Grade III
Contaminants	: Inorganic phosphate ≤0.01mg/g Alkaline phosphatase ≤1 × 10 ⁻⁴ U/mg

PROPERTIES

Stability	: Stable at 4°C for at least 6 months	(Fig.4)
Molecular weight	: approx. 22,000 (SDS-PAGE)	

APPLICATIONS

This peptide fragment is useful as a blocking reagent for ELISA, being an alternative to BSA.

FEATURES

1. Efficiency	: Excellent blocking performance.	(Fig.1)
2. Operation Tim	: Shorter time of operation	(Fig.2)
3. Sensitivity	: No interference on ELISA assay	(Fig.3)

Examples


Fig. 1. BPF-BSA Comparison: Blocking Efficiency

Peroxidase solutions containing BPF-301 or BSA were added to 96 wells polystyrene plate. The plate was incubated at 37°C for 1hr and then washed with 0.02% Tween20. The blocking performance was evaluated by measuring the peroxidase activities retained on polystyrene plate. The lower absorbance shows the higher blocking performance. Blocking performance of BPF-301 is much higher than that of BSA.


Fig. 2. BPF-BSA Comparison: Blocking Time

BPF-301 or BSA solution was added to 96 wells polystyrene plate. The plate containing these blocking reagent was incubated at 4°C for 5~180 minutes, and then washed with 0.02% Tween20. The blocking performance was measured by human serum and anti-human IgG goat antibody POD conjugate system. Operation time required for blocking reaction is shorter than that of BSA.


Fig. 3. BPF-BSA Comparison: Measuring sensitivity

BPF-301 or BSA was used as a blocking reagent in ELISA assay for detection of hCEA. It was shown that no interference on ELISA assay was found in BPF-301.


Fig. 4. Stability (Powder form)

Blocking performance was measured by the method mentioned in Fig. 1.

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