

# dNTP Mixture ( 2mM )

Code No. NTP-201  
Lot No. \*\*\*\*  
Storage Store at -20  
Size 2  $\mu$ moles each

## C e r t i f i c a t e o f A n a l y s

Product : dNTP Mixture ( 2mM )

Appearance : Aqueous solution (Na salt)

pH : 7 . 0

Volume : 1 . 0 ml ( 2 $\mu$ moles each)

Composition	:	Concentratio n	,	Purity(HPLC)
dATP		<u>2 . 0</u>	mM ,	<u>****</u> %
dGTP		<u>2 . 0</u>	mM ,	<u>****</u> %
dCTP		<u>2 . 0</u>	mM ,	<u>****</u> %
dTTP		<u>2 . 0</u>	mM ,	<u>****</u> %